

Limited Warranty

ACCESS INDUSTRIES, 4001 East 138th Street, Grandview, Missouri 64030, warrants to the ORIGINAL PURCHASER ONLY for a **five year period on the drive train** (the drive train is comprised of the motor and gearbox), and a **two year period on all component parts**, [outdoor package is limited to a two year period on the drive train and a one year period on all component parts] commencing with the date of purchase, that each lift and all parts used in the lift shall be free of defects in material and workmanship. The self-addressed Limited Warranty card below must be completed and returned to ACCESS INDUSTRIES within 30 days from the date of the lift's original purchase in order for this Limited Warranty to go into effect. OTHERWISE, THIS LIMITED WARRANTY WILL NOT BE APPLICABLE.

The following are examples of items not considered manufacturing defects and are not covered by this Limited Warranty. 1. Difficulty or damage resulting from improper assembly, operation, care of maintenance. 2. Minor mechanical adjustments such as tightening of nuts, bolts and screws. 3. Damage or difficulties due to misuse, neglect, fire, flood, acts of God, ordinary wear and tear, and damage during shipment that would be the responsibility of the carrier.

NO WARRANTY IS EXTENDED, EXPRESS OR IMPLIED, WHETHER OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, AFTER EXPIRATION OF ONE YEAR FROM THE DATE OF THE ORIGINAL PURCHASE OF THE LIFT. ACCESS INDUSTRIES and its dealer shall not be liable for any consequential, special or incidental damages arising out of the purchase or use of the lift or resulting from the breach of this Limited Warranty, or any implied warranty, and the limit of liability of ACCESS INDUSTRIES and its dealer hereunder shall be the lift's purchase price. Some states do not allow limitations on how long an implied warranty lasts or the exclusion or limitations may not apply to you.

To make a claim under this Limited Warranty, immediately send to the dealer, who sold you the lift and has the responsibility for handling your warranty claim, a notice of your claim, setting forth the defect complained of, along with a copy of the invoice or other proof of purchase and the serial number of the lift. The dealer will notify ACCESS INDUSTRIES. Your notice must be received by the dealer within the warranty period from the date of the lift's original purchase. Any defect in the lift due to materials or workmanship will be repaired or the lift replaced without charge for materials. If the lift is not delivered to the dealer for repair or replacement, the dealer may make a service charge for going to the location of the lift. Such charge is not covered under the Limited Warranty. ALSO, THIS LIMITED WARRANTY DOES NOT COVER ANY LABOR OR SERVICE CHARGES INCURRED IN THE REPAIR OR REPLACEMENT OF A DEFECTIVE LIFT OR PART BY THE DEALER OR ACCESS INDUSTRIES, INC. It is permissible to have any repairs made or replacement work done as a result of any defects in material and workmanship by someone other than the dealer under this Limited Warranty. However, the Limited Warranty does not cover any charges or expenses assessed by any such other person or company performing such repairs or replacement work. All parts used to replace defective materials must be genuine ACCESS INDUSTRIES parts to be covered by this Limited Warranty.

This Limited Warranty gives you specific legal rights, and you may have other rights which vary from state to state.

Learn about new products, services and features at
www.accessind.com



Silver-Glide® HD

Stairway Lift

Owner's Manual & Warranty Information



4001 East 138th Street
Grandview, MO 64030-2837
816-763-3100
www.accessind.com

It is important that this entire manual be read and understood prior to attempting operation of equipment. Also, please observe all cautions and warnings!

Contents

Introduction	3
Operation and Description of Features	4
Keeping your Lift Trouble Free	6
Maintenance	7
Troubleshooting	7
Limited Warranty	8

General Specifications

Rated Load:	375 pounds (170kg)
Rated Speed:	18 feet per minute
Power Supply:	115 VAC, 60 Hz, 3-wire grounded
Motor:	1/2 hp, DC (reversible)
Drive System:	1/8" steel aircraft cable on winding drum
Controls:	Constant pressure switches
Track:	Extruded aluminum
Angle of Incline:	25° to 45°
Application:	Indoor use only

IMPORTANT

The serial number of your Silver-Glide® HD Stairway Lift is:

This serial number contains important information about your lift. This number may be needed if the lift requires replacement parts. Please keep it handy.

Maintenance

Regular maintenance is essential in keeping your Silver-Glide in a proper operating condition. Please remember, you as the owner of this equipment are responsible for making sure maintenance and upkeep is performed on a regularly scheduled basis.

We strongly urge you to protect your investment by having only an authorized Access Industries' dealer perform all required maintenance, service and repair work.

To assure proper operating condition of your lift, the items listed below should be inspected and/or serviced every six (6) months.

- Inspect trailing cable for wear. The cable must be replaced if showing signs of wear. Use talcum powder to lubricate entire length of trailing cable jacket.
- Tighten all fastening anchors.
- Clean the inside of the track to remove dirt and foreign substances. Apply a thin coat of light bodied general purpose grease to the vertical surfaces of the track where the chassis wheels travel. The track should be cleaned and lubricated approximately every 2 months.

Troubleshooting

Before you call for service. . .

A quick check of the Operating Instructions in this manual may prevent an unnecessary service call. Before calling your authorized Access Industries' dealer for service, refer to the following checklist:

If lift will not operate:

- Check for a blown fuse or tripped circuit breaker at the main circuit box.
- Check for obstructions in the path of the lift.
- Check that the swivel seat is in the locked position.
- Replace the 9V batteries in the wireless remote controls if necessary. (make sure the + and - are installed correctly)

If lift will still not operate, contact your local authorized dealer for immediate service. If your lift needs servicing, use only original replacement parts and hardware to maintain quality and performance.

Keeping Your Lift Trouble Free

Your Silver-Glide stairway lift will provide years of trouble free service if you follow these helpful hints:

- Do not exceed capacity of 375 pounds (170kg).
- Do not operate the lift outdoors. The lift is designed for indoor use only.
- Do not operate the lift with a child or pet in your lap. The lift is designed for only one person.
- Do not use an extension cord to supply power to your lift. The lift's power cord should be connected directly to the power source.
- Do not allow weight on lift to be "off center".
- Do not allow children to play on the lift.
- Do not stand up on the lift while operating.
- Do not operate the lift if it has been damaged in any manner.
- Do not attempt to service the lift yourself. Contact your local authorized Access Industries' dealer.
- Do complete and mail in your warranty card or register at www.accessind.com/register
- Do show a family member and/or neighbor how to operate the lift in the event you would need assistance.
- Do keep your lift clean and free of dirt.
- Do have your dealer perform general maintenance and lubrication on your lift to keep it looking and operating like new.
- Do make sure the dealer has explained proper operation and maintenance procedures.

Introduction

Now that you have a new Silver-Glide® HD stairway lift, you can begin enjoying new horizons of independent living!

Your Silver-Glide has been manufactured with the same precision and care that has been a tradition for more than 50 years. You can be proud of the fact that you join thousands of others who have turned to us for the ideal solution to their accessibility problems.

As you will soon discover, your lift is easy to operate and has been designed for years of dependable and trouble free operation. However, as with all electrical and mechanical equipment, proper operation and periodic maintenance is required.

This manual has been prepared to provide clear and concise operating and maintenance instructions. We urge you to read and adhere to the important guidelines provided in this manual before you attempt to operate your lift.

Proper and safe operation of your stairway lift requires proper installation, maintenance and service. Only an authorized Access Industries' dealer should install, maintain, service or repair your new stairway lift. If you need information regarding the local dealer in your area, please call 800-829-9760.

Dealer Information

Name: _____

Address: _____

Telephone No. _____

If you have questions concerning the operation or maintenance of your Silver-Glide, please contact your dealer.

Operation and Description of Features

Learn the location and function of all features before using your Silver-Glide stairway lift. The dealer should explain the features and operating instructions once the lift is installed. If the dealer did not explain this to you, please contact them.

Operating Controls: There is a control switch located underneath one of the armrests on the seat and two wireless remote controls.

To operate the lift, press and hold the switch on the unit control or one of the buttons on a remote control (*red is UP and green is DOWN*).

The lift automatically stops at the ends of the track and it can be stopped at anytime by releasing the switch or button.

Unit Key Lock: An optional key lock limits unwanted use.

Folding Seat: Seat folds up to provide more stairway room when the lift is not in use.

Folding Footrest: The footrest can be folded up and out of the way when the lift is not in use. Also, obstruction sensors are provided in the footrest to detect obstructions. When an obstruction is encountered the lift will stop and the obstruction can be removed. The footrest must be in the down position when riding the lift.

Swivel Seat Handle: Swivel seat handles, provided on each side of the seat, allow the seat to be swiveled 60° or 90° at the upper and lower landings.

Seat Belt: Seat belt must be securely fastened around waist or chest before operating the lift. (not shown)

